


3E'16 Ethan Aw, 2J'16 Lim Jie Yi Joel & 2B'16 Felix Sun Zhiyuan Teacher-Mentor: Mr Ng Kai Yong

2nd Raffles Institution-Singapore Japanese Secondary School Collaboration on Cultural and Science Students' Exchange



An official event of SJ50





Starch Test

- Iodine is used to test for the presence of starch スターチ in food samples
- Iodine ヨウ素 will turn blue-black when added to starch

Which one gives a positive test for starch?






Many plants, such as sweet potato $\mathcal{A}\mathcal{I} - \mathcal{K}$, $\mathcal{F}\mathcal{K}$, yam $\mathcal{K}\mathcal{L}\mathcal{I}$, pumpkin $\mathcal{K}\mathcal{V}\mathcal{I}$, store excess food as starch.

Procedure for Starch Test

1. Add a few drops of iodine solution to a piece of solid food.





Source of photo: https://s3.amazonaws.com/cla ssconnection/951/flashcards/4 208951/jpg/biurets_test-14C0B5226AF61394552.jpg

Biuret Test

Used to test for the presence of proteins タンパク質


Biuret reagent will:

- turn violet in the presence of proteins (+)
- remain blue if there are **no** proteins (-)









Do you know what egg white contains?

Thus, what do you think the colour will be when egg white is tested?

Procedure for Biuret Test

- 1. Add 1 cm³ of sodium hydroxide solution (水酸化ナトリウム) to an equal volume of the food sample solution, mix well.
- 2. Next, add 1% copper (II) sulfate solution (硫酸銅) one drop at a time to the food sample solution.
- 3. Shake well to mix the liquids together.
- 4. Observe and record any colour change.



Benedict's Test

Benedict's reagent is used to test for the presence of reducing sugars 還元糖 in a solution.

Examples of reducing sugars: glucose ブドウ糖, fructose 果糖, galactose ガラクトース, lactose 乳糖, maltose 麦 芽糖

Benedict's reagent changes colour when heated in a mixture with reducing sugar(s).

Procedure for Benedict's Test

- 1. Add 2 cm^3 of the food sample solution to a test tube.
- 2. Add an equal volume of Benedict's solution to the test tube and

swirl the mixture to mix the liquids well.

Benedict's solution

est sampl


water bath

5 minutes and observe a colour change

- 3. Leave the test tube in a <u>boiling</u> water bath for about 5 minutes, or until the colour of the mixture does not change.
- 4. Observe the colour changes during that time as well as the final colour.

With higher concentrations of reducing sugars, the colour of the (heated) mixture will change from light blue to green, yellow, orange, and red, in increasing amount of sugars.

What do you think the colour of the mixture is when Benedict's reagent is <u>heated</u> after being mixed with **apple juice**?



Source of photo: http://brilliantbiologystudent.weebly.com/ethanol-emulsion-test-for-lipids.html

The End Do you have any questions?